

Product datasheet

Anti-IGFBP4 antibody ab4253

2 References

製品の概要

| | |
|----------|---|
| 製品名 | Anti-IGFBP4 antibody |
| 製品の詳細 | Rabbit polyclonal to IGFBP4 |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: IP, WB |
| 種交差性 | 交差種: Mouse, Cow, Human, Pig |
| 免疫原 | HPLC purified IGFBP 4 isolated from conditioned medium from the human glioblastoma cell line, T98G. |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle. |
| 精製度 | Whole antiserum |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab4253** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
| IP | | |
| WB | | |

追加情報

IP: 4µl immunoprecipitated >90% of 10ng IGFBP 4 from 500µl concentrated serum free T98G culture supernatant.
 WB: 1/1000. Detects 1ng SDS-denatured human IGFBP 4 in concentrated serum free T98G culture supernatant. Predicted molecular weight: 24 kDa.
 Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

| | |
|-------|---|
| 機能 | IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors. |
| 配列類似性 | Contains 1 IGFBP N-terminal domain. Contains 1 thyroglobulin type-1 domain. |
| 細胞内局在 | Secreted. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors