

## Product datasheet

# Anti-IL10RA antibody [RM0101-9E12] ab86466

### 製品の概要

製品名	Anti-IL10RA antibody [RM0101-9E12]
製品の詳細	Armenian Hamster monoclonal [RM0101-9E12] to IL10RA
由来種	Armenian hamster
アプリケーション	適用あり: WB
種交差性	交差種: Mouse
免疫原	Recombinant fragment, corresponding to the extracellular domain of Mouse IL10RA

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS
特記事項 (精製)	IgG fraction of culture supernatant purified by Protein A/G affinity chromatography and 0.2 µm filtered.
ポリ/モノ	モノクローナル
クローン名	RM0101-9E12
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab86466** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Predicted molecular weight: 64 kDa.

### ターゲット情報

<b>機能</b>	Receptor for IL10; binds IL10 with a high affinity.
<b>組織特異性</b>	Spleen, thymus, and PBMC. Weak expression in pancreas, skeletal muscle, brain, heart, and kidney. Placenta, lung, and liver showed intermediate levels. Monocytes, B-cells, large granular lymphocytes, and T-cells express high levels.
<b>関連疾患</b>	Inflammatory bowel disease 28
<b>配列類似性</b>	Belongs to the type II cytokine receptor family.
<b>細胞内局在</b>	Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors