

Product datasheet

Anti-IL32 antibody ab62580

2 References 画像数 1

製品の概要

| | |
|--------------|--|
| 製品名 | Anti-IL32 antibody |
| 製品の詳細 | Rabbit polyclonal to IL32 |
| 由来種 | Rabbit |
| 特異性 | ab62580 detects the largest isoform of IL32. |
| アプリケーション | 適用あり: WB |
| 種交差性 | 交差種: Mouse, Human |
| 免疫原 | A 16 amino acid peptide from near the center of human IL32 |
| ポジティブ・コントロール | Jurkat cell lysate. |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | Preservative: 0.02% Sodium Azide Constituents: PBS |
| 精製度 | Immunogen affinity purified |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab62580** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| WB | | Use a concentration of 2.5 - 5 µg/ml. Detects a band of approximately 26 kDa (predicted molecular weight: 26 kDa). |

ターゲット情報

| | |
|-------|---|
| 機能 | Cytokine that may play a role in innate and adaptive immune responses. It induces various cytokines such as TNFA/TNF-alpha and IL8. It activates typical cytokine signal pathways of NF-kappa-B and p38 MAPK. |
| 組織特異性 | Selectively expressed in lymphocytes. Expression more prominent in immune cells than in non-immune cells. |
| 細胞内局在 | Secreted. |

画像



Lane 1 : Anti-IL32 antibody (ab62580) at 2.5 µg/ml

Lane 2 : Anti-IL32 antibody (ab62580) at 5 µg/ml

All lanes : Jurkat cell lysate

Lysates/proteins at 15 µg per lane.

Predicted band size: 26 kDa

Observed band size: 26 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors