

Product datasheet

Anti-IL3RA antibody [6H6] (Allophycocyanin) ab95781

画像数 1

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-IL3RA antibody [6H6] (Allophycocyanin) |
| 製品の詳細 | Mouse monoclonal [6H6] to IL3RA (Allophycocyanin) |
| 標識 | Allophycocyanin. Ex: 645nm, Em: 660nm |
| アプリケーション | 適用あり: Flow Cyt |
| 種交差性 | 交差種: Human |
| 免疫原 | Human IL3RA-transfected COS cells |
| ポジティブ・コントロール | Normal Human peripheral blood cells |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | Preservative: 0.09% Sodium Azide Constituents: PBS, BSA, Gelatin, pH 7.2 |
| 精製度 | Protein G purified |
| ポリ/モノ | モノクローナル |
| クローン名 | 6H6 |
| アイソタイプ | IgG1 |
| 軽鎖の種類 | kappa |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab95781** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------------------------------------|
| Flow Cyt | | Use 5µl for 10 ⁶ cells. |

ターゲット情報

| | |
|-------|---|
| 機能 | This is a receptor for interleukin-3. |
| 配列類似性 | Belongs to the type I cytokine receptor family. Type 5 subfamily. |
| ドメイン | The WSXWS motif appears to be necessary for proper protein folding and thereby efficient intracellular transport and cell-surface receptor binding. The box 1 motif is required for JAK interaction and/or activation. |
| 細胞内局在 | Membrane. |

画像



Flow Cytometry - IL3RA antibody [6H6]
(Allophycocyanin) (ab95781)

Staining of normal human peripheral blood cells with mouse IgG1 k Isotype Control APC (open histogram) or ab95781 (filled histogram). Cells in the lymphocyte gate were used for analysis.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors