


## Product datasheet

# Anti-IL9 antibody ab111915

### 画像数 1

#### 製品の概要

製品名	Anti-IL9 antibody
製品の詳細	Rabbit polyclonal to IL9
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat 
免疫原	Recombinant fragment, corresponding to amino acids 26-144 of Human IL9 (BC066285).
ポジティブ・コントロール	Human fetal liver lysate.

#### 法規制情報

### 医薬用外毒物

#### 製品の特性

製品の状態	Lyophilised: Reconstitute in 200ul Sterile Distilled Water.
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 1% BSA, 98.98% PBS
精製度	Immunogen affinity purified
特記事項(精製)	ab111915 is purified by a peptide affinity column.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab111915** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/2000. Detects a band of approximately 16 kDa (predicted molecular weight: 16 kDa).
IHC-P		1/100 - 1/500.

## ターゲット情報

機能	Supports IL-2 independent and IL-4 independent growth of helper T-cells.
配列類似性	Belongs to the IL-7/IL-9 family.
細胞内局在	Secreted.

## 画像



Anti-IL9 antibody (ab111915) at 1/1000 dilution + Human fetal liver lysate

**Predicted band size:** 16 kDa

**Observed band size:** 16 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors