

Product datasheet

Anti-Infectious Bronchitis (Massachusetts) antibody ab31671

1 References

製品の概要

製品名	Anti-Infectious Bronchitis (Massachusetts) antibody
製品の詳細	Chicken polyclonal to Infectious Bronchitis (Massachusetts)
由来種	Chicken
特異性	Infectious Bronchitis (Massachusetts)
アプリケーション	適用あり: WB, ELISA, Neutralising
免疫原	Infectious Bronchitis (Massachusetts).

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50% for extra stability and aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: Whole serum.
精製度	Whole antiserum
特記事項 (精製)	This antibody is at least 98% pure. No contaminants detected. Single band by SDS PAGE, IEP, and/or RID.
ポリ/モノ	ポリクローナル
アイソタイプ	IgY

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab31671** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution. PubMed: 20534109
ELISA		Use at an assay dependent dilution.

アプリケーション	Abreviews	特記事項
Neutralising		Use at an assay dependent dilution.

ターゲット情報

関連性	Infectious bronchitis is an acute and highly contagious respiratory disease of chickens. The disease is characterized by respiratory signs including gasping, coughing, sneezing, tracheal râles, and nasal discharge. It is caused by a coronavirus. It is an enveloped, single stranded RNA virus. Three virus specific proteins have been identified; the spike glycoprotein, the membrane or matrix glycoprotein, and the nucleocapsid protein.
-----	---

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors