

Product datasheet

Anti-KIR3DL2 antibody - N-terminal ab95303

画像数 1

製品の概要

製品名	Anti-KIR3DL2 antibody - N-terminal
製品の詳細	Rabbit polyclonal to KIR3DL2 - N-terminal
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide selected from the N terminal region of Human KIR3DL2 (NP_006728.2) conjugated to KLH.
ポジティブ・コントロール	HL60 cell line lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
特記事項(精製)	Purified through a protein A column, followed by peptide affinity purification.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab95303** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/100 - 1/500. Predicted molecular weight: 50 kDa.

## ターゲット情報


---

機能	Receptor on natural killer (NK) cells for HLA-A alleles. Inhibits the activity of NK cells thus preventing cell lysis.
配列類似性	Belongs to the immunoglobulin superfamily. Contains 3 Ig-like C2-type (immunoglobulin-like) domains.
細胞内局在	Cell membrane.

---

## 画像

---



Anti-KIR3DL2 antibody - N-terminal (ab95303)  
at 1/100 dilution + HL60 cell line lysate at 35  
µg

**Predicted band size: 50 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors