

Product datasheet

Anti-KRBA2 antibody ab121721

画像数 2

製品の概要

製品名	Anti-KRBA2 antibody
製品の詳細	Rabbit polyclonal to KRBA2
アプリケーション	適用あり: ICC/IF, IHC-P
種交差性	交差種: Human
免疫原	Recombinant Protein Epitope Signature Tag: SVLDSDSGVE FTNQVVHELN ELWPDLKMS GKYHPGQSQG SLEGASRDVK NMISTWMQSN HSCHWAKGLR FMQMVRNQAF DVSL, corresponding to amino acids 320 - 403 of Human KRBA2. Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Human heart muscle RT-4 lysate and Human liver tissue lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab121721** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 1 - 4 µg/ml. Recommend PFA Fixation and Triton X-100 treatment
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報

配列類似性


Contains 1 integrase catalytic domain.

Contains 1 KRAB domain.

翻訳後修飾


Phosphorylated upon DNA damage, probably by ATM or ATR.

画像



Immunofluorescent staining of Human cell line A-431 shows positivity in nuclei but not nucleoli. Recommended concentration of ab121721 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-KRBA2 antibody (ab121721)



ab121721, at a 1/200 dilution, staining KRBA2 in paraffin-embedded Human heart muscle tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-KRBA2 antibody (ab121721)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors