


Product datasheet

Anti-KRT24 antibody ab87195

1 References [画像数 1](#)

製品の概要

製品名	Anti-KRT24 antibody
製品の詳細	Rabbit polyclonal to KRT24
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Horse, Cat, Pig 
免疫原	A synthetic peptide corresponding to a region within internal sequence amino acids 468-517 (SVNMGSRDLVSGDSRSGSCSGQGRDSSKTRVTKTIVEELVDGKVVSSQV S) of Human KRT24.
ポジティブ・コントロール	Human fetal kidney lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab87195** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 55 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

関連性	KRT24 is a member of the type I (acidic) keratin family, which belongs to the superfamily of intermediate filament (IF) proteins. Keratins are heteropolymeric structural proteins which form the intermediate filament. These filaments, along with actin microfilaments and microtubules, comprise the cytoskeleton of epithelial cells. KRT24 is highly expressed in keratinocytes, placenta, colon and spleen, and has been found to be up-regulated in the mucosa of patients suffering of colorectal cancer.
細胞内局在	Intermediate filament; keratin

画像



Anti-KRT24 antibody (ab87195) at 1 µg/ml +
human fetal kidney lysate at 10 µg

Secondary

HRP-conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 55 kDa

Observed band size: 55 kDa

Gel concentration 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors