

Product datasheet


Kainic acid ab120100

39 References

製品の概要

<b>製品名</b>	Kainic acid
<b>製品の詳細</b>	Prototypic kainate receptor agonist
<b>生理活性の詳細</b>	Prototypic agonist at the kainate class of ionotropic glutamate receptors. Potent excitant and neurotoxin, used to model epilepsy and neurodegenerative states. Also available in Kit(s): Ionotropic agonists ( <a href="#">ab120323</a> ).
<b>精製度</b>	Also available in simple stock solutions ( <a href="#">ab144490</a> ) - add 1 ml of water to get an exact, ready-to-use concentration. > 99%

製品の特性

<b>体系名</b>	(2S,3S,4S)-3-(Carboxymethyl)-4-(prop-1-en-2-yl)pyrrolidine-2-carboxylic acid
<b>分子量</b>	213.23
<b>構造式</b>	
<b>分子式</b>	C <sub>10</sub> H <sub>15</sub> NO <sub>4</sub>
<b>CAS 番号</b>	487-79-6
<b>保存方法</b>	Store at +4°C. Store under desiccating conditions. The product can be stored for up to 12 months.
<b>溶解性</b>	Soluble in water to 25 mM
<b>使用に関する注意</b>	Wherever possible, you should prepare and use solutions on the same day. However, if you need to make up stock solutions in advance, we recommend that you store the solution as aliquots in tightly sealed vials at -20°C. Generally, these will be useable for up to one month. Before use, and prior to opening the vial we recommend that you allow your product to equilibrate to room temperature for at least 1 hour.  Need more advice on solubility, usage and handling? Please visit our <a href="#">frequently asked questions (FAQ) page</a> for more details.
<b>SMILES 線形表記</b>	CC(=C)[C@H]1CNC([C@H]1CC(=O)O)C(=O)O
<b>由来</b>	<i>Digenea simplex</i>

## アプリケーション

---

Our [Abpromise guarantee](#) covers the use of **ab120100** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Functional Studies		Use at an assay dependent concentration.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE, NOT FOR USE IN HUMANS"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors
- Abcam biochemicals are novel compounds and we have not tested their biological activity in house. Please use the literature to identify how to use these products effectively. If you require further assistance please contact the scientific support team