

## Product datasheet

# Anti-LRRN5 antibody ab125921

1 References [画像数 2](#)

### 製品の概要

---

製品名	Anti-LRRN5 antibody
製品の詳細	Rabbit polyclonal to LRRN5
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB, IHC-P
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment, corresponding to a region within amino acids 213-535 of Human LRRN5 (O75325).
ポジティブ・コントロール	U-87 MG whole cell lysate; Human Hepatoma tissue.

### 法規制情報

---

## 医薬用外毒物

### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.75% Glycine, 1.21% Tris, 20% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

---

Our [Abpromise guarantee](#) covers the use of **ab125921** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 79 kDa.
IHC-P		1/100 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## ターゲット情報

**関連性** LRRN5 (Leucine rich repeat neuronal 5) is a transmembrane protein belonging to the leucine-rich repeat (LRR) superfamily. It has homology with other proteins that function as cell-adhesion molecules or as signal transduction receptors. It is overexpressed in malignant gliomas.

**細胞内局在** Membrane; Single-pass type I membrane protein.


## 画像



Anti-LRRN5 antibody (ab125921) at 1/1000 dilution + U-87 MG whole cell lysate at 30 µg

**Predicted band size:** 79 kDa

Gel concentration: 7.5%



ab125921, at 1/500, staining LRRN5 in Paraffin Embedded Human Hepatoma tissue by Immunohistochemistry.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors