


Product datasheet

Anti-MAL2 antibody ab75347

2 References

製品の概要

製品名	Anti-MAL2 antibody
製品の詳細	Rabbit polyclonal to MAL2
アプリケーション	適用あり: ELISA, ICC, WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat 
免疫原	Synthetic peptide derived from the N-terminal domain of human MAL2

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab75347** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
ICC		1/50 - 1/200.
WB		1/200 - 1/1000. Predicted molecular weight: 19 kDa.

ターゲット情報

機能	Member of the machinery of polarized transport. Required for the indirect transcytotic route at the step of the egress of the transcytosing cargo from perinuclear endosomes in order for it to travel to the apical surface via a raft-dependent pathway.
組織特異性	Predominantly expressed in kidney, lung, and liver. Also found in thyroid gland, stomach and, at lower levels in testis and small intestine.
配列類似性	Belongs to the MAL family. Contains 1 MARVEL domain.
細胞内局在	Cell membrane. Apical cell membrane. Endomembrane system. Cytoplasm > perinuclear region. Associated with lipid rafts. In polarized epithelial cells, restricted to the apical surface. In hepatocytes, as well as in polarized hepatoma HepG2 cells, found in the canalicular membrane, equivalent to the apical surface, beneath the canalicular actin cytoskeleton. In non-polarized HepG2 cells, distributed to the perinuclear region.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors