

## Product datasheet

# Mouse monoclonal [MARE-1] Anti-Rat IgE heavy chain (FITC) ab11667

### 製品の概要

製品名	Mouse monoclonal [MARE-1] Anti-Rat IgE heavy chain (FITC)
由来種	Mouse
ターゲット生物種	Rat
特異性	This antibody is specific for the epsilon chain of rat immunoglobulin IgE and does not cross-react with other classes of rat immunoglobulin.
アプリケーション	<b>適用あり:</b> Flow Cyt
免疫原	Rat IgE myeloma proteins (IR162, IR1016, IR2 and IR410).
標識	FITC. Ex: 493nm, Em: 528nm

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 50% Glycerol, PBS, pH 7.4
精製度	Immunogen affinity purified
特記事項 (精製)	Purified IgG prepared by affinity chromatography.
ポリ/モノ	モノクローナル
クローン名	MARE-1
ミエローマ	Sp2/0
アイソタイプ	IgG1
特記事項	This product should be stored undiluted. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab11667** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

Flow Cyt

Use a concentration of 5 µg/ml.

Use 50 µl of the suggested working dilution to label 2 x 10<sup>7</sup> cells in 100 µl.

[ab91356](#) - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors