

Product datasheet

Mouse monoclonal [MARG1-2] Anti-Rat IgG1 heavy chain (Biotin) ab11668

★★★★★ 1 Abreviews 1 References

製品の概要

製品名	Mouse monoclonal [MARG1-2] Anti-Rat IgG1 heavy chain (Biotin)
由来種	Mouse
ターゲット生物種	Rat
特異性	Ab11668 recognises rat IgG1. No reactivity with rat IgG2a, IgG2b, IgG2c, IgA, IgM, IgE or IgD.
アプリケーション	適用あり: ELISA
免疫原	Rat IR27 IgG1 myeloma protein.
標識	Biotin

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 50% Glycerol, PBS, pH 7.4
特記事項 (精製)	Purified IgG prepared by affinity chromatography.
ポリ/モノ	モノクローナル
クローン名	MARG1-2
ミエローマ	Sp2/0
アイソタイプ	IgG1
特記事項	

Immunoglobulins, or antibodies, are a complex, heterogeneous mixture of proteins that exhibit two fundamental types of structural variation. Subtle structural differences in their antigen combining sites, or variable regions, account for their unique antigen binding specificities.

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab11668** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use a concentration of 0.5 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors