

## Product datasheet

# Anti-NCAM antibody [B-A19] ab27299

1 References [画像数 1](#)

### 製品の概要

製品名	Anti-NCAM antibody [B-A19]
製品の詳細	Mouse monoclonal [B-A19] to NCAM
アプリケーション	<b>適用あり:</b> Flow Cyt
種交差性	<b>交差種:</b> Human
免疫原	Tissue/ cell preparation of U937 cell line (Human)
ポジティブ・コントロール	human lymphocytes

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium Azide Constituents: 1% BSA, PBS
精製度	Ion Exchange Chromatography
ポリ/モノ	モノクローナル
クローン名	B-A19
ミエローマ	x63-Ag8.653
アイソタイプ	IgG1

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab27299** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

Flow Cyt

### 追加情報

Flow Cyt: Use 10µl for 10<sup>6</sup> cells or 100µl of whole blood.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報


---

機能	This protein is a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc.
配列類似性	Contains 2 fibronectin type-III domains. Contains 5 Ig-like C2-type (immunoglobulin-like) domains.
細胞内局在	Secreted and Cell membrane.

---

## 画像

---



A typical staining pattern with the B-A19 monoclonal antibody of lymphocytes FACS analysis of NCAM expression in human lymphocytes using ab27299.

Flow Cytometry / FACS - NCAM antibody [B-A19]  
(ab27299)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors