

Product datasheet

Anti-NCAM1 antibody [B-A19] ab27299

1 References [画像数 1](#)

製品の概要

製品名	Anti-NCAM1 antibody [B-A19]
製品の詳細	Mouse monoclonal [B-A19] to NCAM1
由来種	Mouse
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Human
免疫原	Tissue/ cell preparation of U937 cell line (Human)
ポジティブ・コントロール	human lymphocytes

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium Azide Constituents: 1% BSA, PBS
精製度	Ion Exchange Chromatography
ポリ/モノ	モノクローナル
クローン名	B-A19
ミエローマ	x63-Ag8.653
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab27299** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

Flow Cyt

追加情報

Flow Cyt: Use 10µl for 10⁶ cells or 100µl of whole blood.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

機能	This protein is a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc.
配列類似性	Contains 2 fibronectin type-III domains. Contains 5 Ig-like C2-type (immunoglobulin-like) domains.
細胞内局在	Secreted and Cell membrane.

画像



A typical staining pattern with the B-A19 monoclonal antibody of lymphocytes FACS analysis of NCAM expression in human lymphocytes using ab27299.

Flow Cytometry - Anti-NCAM1 antibody [B-A19]
(ab27299)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors