

## Product datasheet

# Anti-Oligodendrocyte myelin glycoprotein antibody ab96426

★★★★☆ 1 Abreviews [画像数 1](#)

### 製品の概要

製品名	Anti-Oligodendrocyte myelin glycoprotein antibody
製品の詳細	Rabbit polyclonal to Oligodendrocyte myelin glycoprotein
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Cow, Macaque monkey, Orangutan
免疫原	Synthetic peptide corresponding to a region within amino acids 249 - 311 of Human Oligodendrocyte myelin glycoprotein.
ポジティブ・コントロール	A549, IMR32 or HCT116 cell lysate.

### 法規制情報

**医薬用外毒物**

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 1 x PBS, pH 7.0
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab96426** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB	★★★★☆	1/500 - 1/3000. Predicted molecular weight: 49 kDa.

## ターゲット情報

機能	Cell adhesion molecule contributing to the interactive process required for myelination in the central nervous system.
組織特異性	Oligodendrocytes and myelin of the central nervous system.
配列類似性	Contains 8 LRR (leucine-rich) repeats. Contains 1 LRRNT domain.
翻訳後修飾	O-glycosylated in its Ser/Thr-rich repeat domain.
細胞内局在	Cell membrane.

## 画像



Anti-Oligodendrocyte myelin glycoprotein antibody (ab96426) at 1/1000 dilution + IMR32 whole cell lysate at 30 µg

**Predicted band size:** 49 kDa

10% SDS PAGE

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors