

Product datasheet

Anti-PBF1 antibody ab98964

画像数 1

製品の概要

製品名	Anti-PBF1 antibody
製品の詳細	Rabbit polyclonal to PBF1
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Arabidopsis thaliana
免疫原	Full length protein produced in E.coli.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab98964** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/1000. Predicted molecular weight: 25 kDa.

ターゲット情報


関連性 The proteasome is a multicatalytic proteinase complex which is characterized by its ability to cleave peptides with Arg, Phe, Tyr, Leu, and Glu adjacent to the leaving group at neutral or

slightly basic pH. The proteasome has an ATP-dependent proteolytic activity.

細胞内局在

Cytoplasmic and Nuclear

画像



Anti-PBF1 antibody (ab98964) at 1/1000 dilution + Arabidopsis crude plant extracts.

Predicted band size: 25 kDa

Observed band size: 25 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors