

Product datasheet

Anti-PCDHA9 antibody ab110512

画像数 1

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-PCDHA9 antibody |
| 製品の詳細 | Rabbit polyclonal to PCDHA9 |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: WB |
| 種交差性 | 交差種: Human |
| 免疫原 | Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 233-261 of Human PCDHA9 (NP_054724.1, NP_114063.1) |
| ポジティブ・コントロール | NCI-H460 cell line lysates |

製品の特性

| | |
|----------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term. |
| バッファー | Preservative: 0.09% Sodium azide Constituent: PBS |
| 精製度 | Immunogen affinity purified |
| 特記事項(精製) | This antibody is purified through a protein A column, followed by peptide affinity purification. |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab110512** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| WB | | 1/100 - 1/500. Predicted molecular weight: 102 kDa. |

| | |
|-------|--|
| 関連性 | Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain. |
| 細胞内局在 | Cell membrane; Single-pass type I membrane protein |

画像



Anti-PCDHA9 antibody (ab110512) at 1/100 dilution + NCI-H460 cell line lysates at 35 µg

Predicted band size: 102 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors