

Product datasheet

Anti-PDE4C antibody ab14608

★★★★☆ 1 Abreviews

製品の概要

製品名	Anti-PDE4C antibody
製品の詳細	Rabbit polyclonal to PDE4C
特異性	This antibody detects 3 isoforms of PDE4C subclasses (PDE4C1, PDE4C2 and PDE4C3).
アプリケーション	適用あり: IP, WB
種交差性	交差種: Human
免疫原	Synthetic cyclic peptide (carboxy-terminal peptide common to all PDE4C proteins).
ポジティブ・コントロール	recombinant protein PDE4C1/2

法規制情報

医薬用外毒物

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Thimerosal / Merthiolate
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab14608** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

アプリケーション	Abreviews	特記事項
IP		1/20.
WB	★★★★☆	1/500. Predicted molecular weight: 80 kDa. Can be blocked with PDE4C peptide (ab169990) .

ターゲット情報

機能	Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes.
組織特異性	Expressed in various tissues but not in cells of the immune system.
パスウェイ	Purine metabolism; 3',5'-cyclic AMP degradation; AMP from 3',5'-cyclic AMP: step 1/1.
配列類似性	Belongs to the cyclic nucleotide phosphodiesterase family. PDE4 subfamily.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors