

Product datasheet

Anti-PF4V1 antibody ab97799

画像数 1

製品の概要

製品名	Anti-PF4V1 antibody
製品の詳細	Rabbit polyclonal to PF4V1
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to a region within amino acids 40-104 (GDLQCLCVKT TSQVRPRHIT SLEVIKAGPH CPTAQLIATL KNGRKICLDL QALLYKKIIK EHLES) of Human PF4V1 (NP_002611) <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	A431, H1299, HeLa or HepG2 lysates

法規制情報

**医薬用外毒物**

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
精製度	Immunogen affinity purified
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab97799** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 12 kDa.

## ターゲット情報

関連性	PF4V1 is an inhibitor of angiogenesis and an inhibitor of endothelial cell chemotaxis (in vitro).
細胞内局在	Secreted

## Anti-PF4V1 antibody 画像



Anti-PF4V1 antibody (ab97799) at 1/1000 dilution + A431 whole cell lysate at 30 µg

**Predicted band size : 12 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors