

Product datasheet

Anti-POLD3 antibody ab74307

画像数 1

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-POLD3 antibody |
| 製品の詳細 | Rabbit polyclonal to POLD3 |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: WB 適用なし: IP |
| 種交差性 | 交差種: Human 交差が予測される動物種: Rabbit, Dog, Pig, Chimpanzee, Rhesus monkey, Gorilla, Orangutan, Bat |
| 免疫原 | Synthetic peptide corresponding to a region between residues 416 and 466 of human POLD3 (NP_006582.1) |
| ポジティブ・コントロール | Lysate from HeLa or 293T cells |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| バッファー | Preservative: 0.09% Sodium Azide Constituents: 8mM PBS, 60mM Citrate, 150mM Tris, pH 7-8 |
| 精製度 | Immunogen affinity purified |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab74307** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| WB | | 1/2000 - 1/10000. Predicted molecular weight: 51 kDa. |

追加情報

Is unsuitable for IP.

ターゲット情報


機能

Required for optimal DNA polymerase delta activity.

細胞内局在

Nucleus.

画像



All lanes : Anti-POLD3 antibody (ab74307) at 0.1 µg/ml

Lane 1 : Whole cell lysate from HeLa cells at 50 µg

Lane 2 : Whole cell lysate from HeLa cells at 15 µg

Lane 3 : Whole cell lysate from HeLa cells at 5 µg

Lane 4 : Whole cell lysate from 293T cells at 50 µg

Predicted band size: 51 kDa

Observed band size: 63 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors