

Product datasheet

Anti-PSTPIP2 antibody ab76793

画像数 2

製品の概要

製品名	Anti-PSTPIP2 antibody
製品の詳細	Mouse polyclonal to PSTPIP2
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein, corresponding to amino acids 1-312 of Human PSTPIP2
ポジティブ・コントロール	transfected 293T cell line, human pancreas

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS, pH 7.2
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab76793** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml. Detects a band of approximately 39 kDa (predicted molecular weight: 39 kDa).

ターゲット情報

関連性 PSTPIP2 (Proline-serine-threonine phosphatase-interacting protein 2) binds to F-actin and may be involved in regulation of the actin cytoskeleton.

細胞内局在 Cell Membrane; peripheral membrane protein and Cytoplasmic (by similarity).

画像




Anti-PSTPIP2 antibody (ab76793) at 1/500 dilution + human pancreas at 25 µg

Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size: 39 kDa

Observed band size: 39 kDa



All lanes : Anti-PSTPIP2 antibody (ab76793) at 1/500 dilution

Lane 1 : PSTPIP2 transfected 293T cell line

Lane 2 : non transfected 293T cell line

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size: 39 kDa

Observed band size: 30 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors