


Product datasheet

Anti-PURA antibody ab125200

1 Abreviews [画像数 1](#)

製品の概要

製品名	Anti-PURA antibody
製品の詳細	Rabbit polyclonal to PURA
由来種	Rabbit
アプリケーション	適用あり: IP 適用なし: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Sheep, Rabbit, Guinea pig, Cow, Dog, Pig, Chimpanzee, Ferret, Rhesus monkey, Gorilla, Orangutan 
免疫原	Synthetic peptide, corresponding to a region within amino acids 100-150 of Human PURA (NP_005850.1).
ポジティブ・コントロール	HeLa whole cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituent: 99% Tris citrate/phosphate Note: pH 7 to 8.
精製度	Immunogen affinity purified
特記事項(精製)	ab125200 was affinity purified using an epitope specific to PURA immobilized on a solid support.
ポリモノ	ポリクローナル
アイソタイプ	IgG


アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab125200** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IP		Use at 2-10 µg/mg of lysate.
追加情報	Is unsuitable for WB.	
ターゲット情報		
機能	This is a probable transcription activator that specifically binds the purine-rich single strand of the PUR element located upstream of the MYC gene. May play a role in the initiation of DNA replication and in recombination.	
配列類似性	Belongs to the PUR DNA-binding protein family.	
細胞内局在	Nucleus.	

画像



Detection of PURA by Western Blot of Immunoprecipitate. Lane 1; an antibody to PURA at 1 µg/ml staining PURA in HeLa whole cell lysate immunoprecipitated using ab125200 at 6 µg/mg lysate (1 mg/IP; 20% of IP loaded/lane). The lysate in lane 2 is precipitated with an IgG control antibody. Detection: Chemiluminescence with exposure time of 10 seconds.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors