

Product datasheet

Anti-Peripherin antibody ab53051

★★★★★ 1 Abreviews [画像数 2](#)

製品の概要

製品名	Anti-Peripherin antibody
製品の詳細	Rabbit polyclonal to Peripherin
由来種	Rabbit
特異性	ab53051 detects endogenous levels of total Peripherin protein.
アプリケーション	適用あり: WB, IHC-P, ELISA
種交差性	交差種: Human 交差が予測される動物種: Mouse
免疫原	A synthetic peptide derived from human Peripherin.
ポジティブ・コントロール	Human brain tissue and HepG2 cell extracts.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4
精製度	Immunogen affinity purified
特記事項 (精製)	ab53051 was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab53051** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 54 kDa (predicted molecular weight: 54 kDa).

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★	Use at an assay dependent concentration.
ELISA		1/5000.

ターゲット情報

機能	Class-III neuronal intermediate filament protein.
配列類似性	Belongs to the intermediate filament family.

画像



Western blot - Peripherin antibody (ab53051)


All lanes : Anti-Peripherin antibody (ab53051) at 1/500 dilution

Lane 1 : HepG2 cell extract

Lane 2 : HepG2 cell extract with immunising peptide

Predicted band size: 54 kDa

Observed band size: 54 kDa



Immunohistochemistry (Paraffin-embedded sections)
- Peripherin antibody (ab53051)

ab53051 at 1/50 dilution staining Peripherin in human brain by Immunohistochemistry, Paraffin embedded tissue (in the absence and presence of the immunising peptide).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors