

## Product datasheet

# Anti-Pleckstrin antibody ab17020

4 References [画像数 1](#)

### 製品の概要

製品名	Anti-Pleckstrin antibody
製品の詳細	Goat polyclonal to Pleckstrin
アプリケーション	<b>適用あり:</b> ELISA, WB
種交差性	<b>交差種:</b> Human <b>交差が予測される動物種:</b> Mouse, Dog
免疫原	Synthetic peptide: EWIKAIQMASRTGK, corresponding to C terminal amino acids 337-350 of Human Pleckstrin. <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	Human heart and PBMC lysates

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 0.5% BSA, Tris buffered saline. pH 7.3
精製度	Immunogen affinity purified
特記事項 (精製)	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
ポリモノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab17020** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		

WB

## 追加情報

ELISA: 1/16000.

WB: Use at a concentration of 0.05 - 1 µg/ml. Detects a band of approximately 45 kDa (predicted molecular weight: 40 kDa) in Human Heart and PBMC lysates

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

## 機能


Major protein kinase C substrate of platelets.

## 配列類似性

Contains 1 DEP domain.

Contains 2 PH domains.

## 画像



Predicted band size : 40 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors