

## Product datasheet

# Anti-Pro + Mature Surfactant Protein B antibody ab40876

★★★★★ 2 Abreviews 8 References 画像数 2

### 製品の概要

製品名	Anti-Pro + Mature Surfactant Protein B antibody
製品の詳細	Rabbit polyclonal to Pro + Mature Surfactant Protein B
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA, IHC-P
種交差性	交差種: Mouse, Sheep, Cow, Human
免疫原	Full length native SP-B protein (purified) (Human)

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium azide Constituent: Whole serum
精製度	Whole antiserum
ポリモノ	ポリクローナル
アイソタイプ	IgG


### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab40876** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB	★★★★★	1/5000. Predicted molecular weight: 42 kDa. (reduced or non-reduced samples).
ELISA		1/500 - 1/1000.
IHC-P		Use at an assay dependent concentration.

画像



Adult mouse lung was embedded in paraffin and stained with ab40876 (1/500 dilution).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Pro + Mature Surfactant Protein B antibody (ab40876)



**All lanes :** Anti-Pro + Mature Surfactant Protein B antibody (ab40876) at 1/5000 dilution

**Lane 1 :** Adult mouse lung homogenate

**Lane 2 :** Human lung homogenate

**Lane 3 :** Surfactant Protein B peptide

**Lane 4 :** HEK-293T lysate

**Predicted band size:** 42 kDa

Western blot - Anti-Pro + Mature Surfactant Protein B antibody (ab40876)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors