

Product datasheet

Anti-Pygopus 2 antibody [PYGAD53A] ab51367

画像数 1

製品の概要

製品名	Anti-Pygopus 2 antibody [PYGAD53A]
製品の詳細	Mouse monoclonal [PYGAD53A] to Pygopus 2
由来種	Mouse
アプリケーション	適用あり: Dot blot, ELISA, WB
種交差性	交差種: Human
免疫原	Recombinant fragment (Human)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
精製度	Protein G purified
特記事項 (精製)	Filtered through a 0.22 µm membrane.
ポリモノ	モノクローナル
クローン名	PYGAD53A
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab51367** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Dot blot		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.


アプリケーション	Abreviews	特記事項
----------	-----------	------

WB Use at an assay dependent dilution. Predicted molecular weight: 41 kDa.

ターゲット情報

機能	Involved in signal transduction through the Wnt pathway.
配列類似性	Contains 1 PHD-type zinc finger.
細胞内局在	Nucleus.

画像



Lane 1 : Molecular weight markers

Lane 2 : Anti-Pygopus 2 antibody
[PYGAD53A] (ab51367)

Lane 2 : immunizing peptide

Predicted band size: 41 kDa

Observed band size: 35,38 kDa

Western blot - Anti-Pygopus 2 antibody
[PYGAD53A] (ab51367)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors