

Product datasheet

Anti-RASA3 antibody ab31286

製品の概要

製品名 Anti-RASA3 antibody
 製品の詳細 Goat polyclonal to RASA3

i This product is a **fast track antibody**. It has been affinity purified and shows high titre values against the immunizing peptide by ELISA.

[Read the terms of use »](#)

種交差性

交差が予測される動物種: Mouse, Rat, Dog, Human 

免疫原

Synthetic peptide:

C-QYKRDKFKKTKYGSQ

, corresponding to C terminal amino acids 796-810 of Human RASA3

 [Run BLAST with](#)

 [Run BLAST with](#)

製品の特性

製品の状態 Liquid
 保存方法 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
 バッファー Preservative: 0.02% Sodium Azide
 Constituents: 0.5% BSA, Tris buffered saline. pH 7.3
 精製度 Immunogen affinity purified
 特記事項 (精製) Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
 ポリ/モノ ポリクローナル
 アイソタイプ IgG

アプリケーション

追加情報 This antibody gave a positive result in ELISA against the immunizing peptide. Antibody detection limit dilution 1:32,000.
 Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

機能	Inhibitory regulator of the Ras-cyclic AMP pathway. Binds inositol tetrakisphosphate (IP4) with high affinity. Might be a specific IP4 receptor.
配列類似性	Contains 1 Btk-type zinc finger. Contains 2 C2 domains. Contains 1 PH domain. Contains 1 Ras-GAP domain.
細胞内局在	Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors