

## Product datasheet

# Rabbit Anti-Cow IgA H&L (Alkaline Phosphatase) ab112627

### 製品の概要

製品名	Rabbit Anti-Cow IgA H&L (Alkaline Phosphatase)
製品の詳細	Rabbit polyclonal Secondary Antibody to Cow IgA (AP)
由来種	Rabbit
ターゲット生物種	Cow
アプリケーション	<b>適用あり:</b> WB, ELISA, IHC-P, ICC/IF
標識	Alkaline Phosphatase

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.10 Preservative: 0.09% Sodium azide Constituents: 1.19% HEPES, 0.58% Sodium chloride, 0.01% Magnesium chloride, 0.001% Zinc chloride, 0.2% BSA
精製度	Immunogen affinity purified
特記事項 (精製)	Antiserum was solid phase adsorbed to ensure class specificity. ab112627 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to alkaline phosphatase.
ポリモノ	ポリクローナル
アイソタイプ	IgA

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab112627** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/10000. chemiluminescent detection.
ELISA		1/1000 - 1/10000. as a primary.

アプリケーション	Abreviews	特記事項
IHC-P		1/100 - 1/1000.
ICC/IF		1/100 - 1/1000.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors