

Product datasheet

Rabbit Anti-Goat IgG H&L (Alkaline Phosphatase) ab6742

2 References 画像数 1

製品の概要

製品名	Rabbit Anti-Goat IgG H&L (Alkaline Phosphatase)
由来種	Rabbit
ターゲット生物種	Goat
アプリケーション	適用あり: Dot blot, ELISA, ICC/IF, WB, IHC-P, IHC-Fr
免疫原	Goat IgG whole molecule
標識	Alkaline Phosphatase

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 8.00 Preservative: 0.1% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.788% Tris HCl, 0.0095% Magnesium chloride, 0.0014% Zinc chloride, 50% Glycerol
精製度	Affinity purified
特記事項 (精製)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Goat IgG coupled to agarose beads.
特記事項 (標識)	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab6742** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Dot blot		Use at an assay dependent concentration.

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.
ICC/IF		1/200 - 1/1000.
WB		Use at an assay dependent dilution.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

画像



Ab6742 was used at dilution 1/4000 with the primary antibody [ab20773](#) in ELISA. See the review on [ab20773](#).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors