

Product datasheet

Rabbit Anti-Pig IgG H&L (FITC) ab6773

1 References

製品の概要

製品名	Rabbit Anti-Pig IgG H&L (FITC)
由来種	Rabbit
ターゲット生物種	Pig
アプリケーション	適用あり: Flow Cyt, ICC/IF, ELISA
免疫原	pig whole molecule
標識	FITC. Ex: 493nm, Em: 528nm

法規制情報

医薬用外毒物

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.60 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.164% Sodium phosphate, 0.878% Sodium chloride, 1% BSA
精製度	Affinity purified
特記事項 (精製)	Diaminotriazinylaminofluorescein (DTAF) (Molecular Weight 530 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 2.3 moles DTAF per mole of Rabbit IgG
特記事項 (標識)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using pig IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
ポリ/モノ	ポリクローナル

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab6773** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		1/500 - 1/2500.
ICC/IF		1/1000 - 1/5000.
ELISA		1/10000 - 1/50000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors