

Product datasheet

Rabbit monoclonal [M1gG51-4] Anti-Mouse IgG1 H&L ab125913

リコンビナント RabMAb[®]

5 References 画像数 1

製品の概要

製品名	Rabbit monoclonal [M1gG51-4] Anti-Mouse IgG1 H&L
製品の詳細	Rabbit monoclonal [M1gG51-4] Anti-Mouse IgG1 H&L
由来種	Rabbit
ターゲット生物種	Mouse
アプリケーション	適用あり: ELISA
免疫原	Synthetic peptide within Mouse IgG1. The exact sequence is proprietary.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.005% Sodium azide Constituents: 99% PBS, 0.1% BSA
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	M1gG51-4
アイソタイプ	IgG
特記事項	Human, Rat: We have preliminary internal testing data to indicate this antibody may not react with these species. Please contact us for more information. Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents This product is a recombinant rabbit monoclonal antibody.

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab125913** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

ELISA

1/1000.

画像



Direct ELISA antibody dose-response curve using ab125913 between 0-1000 ng/ml.

Antigen concentration of 1000 ng/mL. An alkaline phosphatase-conjugated goat anti-rabbit IgG (H+L) (1/2500) was used as the secondary antibody.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors