

Product datasheet

Rat IgG2b, kappa monoclonal [RTK4530] (PE/Cy5®) - Isotype control ab18538

1 References

製品の概要

| | |
|-----------|--|
| 製品名 | Rat IgG2b, kappa monoclonal [RTK4530] (PE/Cy5®) - Isotype control |
| 標識 | PE/Cy5®. Ex: 496nm, Em: 670nm |
| 特記事項 (標識) | The solution is free of unconjugated PE/Cy5 and unconjugated antibody. |
| アプリケーション | 適用あり: Flow Cyt |
| 免疫原 | Trinitrophenol and KLH |
| 特記事項 | Isotype control has been screened on a variety of resting, activated, live, and fixed mouse and human tissues. |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS |
| 精製度 | Immunogen affinity purified |
| ポリ/モノ | モノクローナル |
| クローン名 | RTK4530 |
| アイソタイプ | IgG2b |
| 軽鎖の種類 | kappa |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab18538** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
| Flow Cyt | | |

追加情報

FACS: Use at an assay dependent dilution for surface or intracellular staining.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors