

Product datasheet

Rat monoclonal [LO-MM-9] Anti-Mouse IgM heavy chain (Biotin) ab11659

製品の概要

製品名	Rat monoclonal [LO-MM-9] Anti-Mouse IgM heavy chain (Biotin)
由来種	Rat
ターゲット生物種	Mouse
特異性	ab11659 recognises the mu chain of mouse IgM.
アプリケーション	適用あり: ELISA
免疫原	Polyclonal and monoclonal mouse IgM.
標識	Biotin

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 50% Glycerol, PBS, pH 7.2
精製度	Purified IgM
特記事項 (精製)	Purified IgG prepared by affinity chromatography on mouse IgM/Sepharose 4B.
ポリ/モノ	モノクローナル
クローン名	LO-MM-9
ミエローマ	IR983F
アイソタイプ	IgG2a
特記事項	This product should be stored undiluted. Storage in frost free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use.

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab11659** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション

Abreviews

特記事項

ELISA

Use a concentration of 0.25 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors