

Product datasheet

Recombinant Cow IP10 protein ab87463

製品の概要

製品名	Recombinant Cow IP10 protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Yeast

アミノ酸配列

生物種	Cow
配列	VPLSRNTRCS CIEISNGSVN PRSLEKLEVI PASQSCPRVE IIATMKKNGE KRCLNPESKT IKNLLKAINK QRTKRSPRTR KEA
領域	20 to 102

特性

Our [Abpromise guarantee](#) covers the use of **ab87463** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilised

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Preservative: None Constituents: 10% Trehalose, 1.0X PBS
再構成	Reconstitute with sterile phosphate-buffered saline containing at least 0.1% carrier protein.

関連情報

機能	Chemotactic for monocytes and T-lymphocytes. Binds to CXCR3.
----	--

配列類似性	Belongs to the intercrine alpha (chemokine CxC) family.
翻訳後修飾	CXCL10(1-73) is produced by proteolytic cleavage after secretion from keratinocytes.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors