

Product datasheet

Recombinant HSV2 gG Envelope protein ab43047

製品の概要

製品名	Recombinant HSV2 gG Envelope protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Escherichia coli

アミノ酸配列

配列の追加情報	The protein contains the HSV-2 gG immunodominant regions.
---------	---

特性

Our [Abpromise guarantee](#) covers the use of **ab43047** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

生理活性	Immunoreactive with sera of HSV-infected individuals, excellent antigen for detection of HSV with minimal specificity problems.
アプリケーション	ELISA Western blot
精製度	> 95 % SDS-PAGE. Sephacrose derived purification.
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Preservative: None Constituents: 50% Glycerol, 25mM Tris HCl, 1mM EDTA
------------	--

関連情報

関連性	HSV2 is a member of the herpes virus family, Herpesviridae. All herpes viruses are
-----	--

morphologically identical. They have a large double stranded DNA genome. The virion consists of an icosahedral nucleocapsid which is surrounded by a lipid bilayer envelope. Between the capsid and the envelope is an amorphous layer of proteins, termed the tegument. Glycoprotein G (gG) is one of seven envelope glycoproteins in HSV2.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors