

Product datasheet

Recombinant Hongotoxin1 protein ab57647

製品の概要

製品名	Recombinant Hongotoxin1 protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli
アミノ酸配列	
生物種	Centruroides limbatus
配列	TVIDVKCTSPKQCLPPCKAQFGIRAGAKCMNGKCKCYPH

特性

Our [Abpromise guarantee](#) covers the use of **ab57647** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

製品の状態	Lyophilised
-------	-------------

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
再構成	Soluble in aqueous buffer.

関連情報

関連性	Hongotoxin1 is a 39aa peptide recently isolated from the venom of <i>Centruroides limbatus</i> , which blocks the voltage-gated K ⁺ channels Kv1.1, Kv1.2, and Kv1.3 at picomolar toxin concentrations.
-----	--

細胞内局在	Secreted
-------	----------

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors