

## Product datasheet

# Recombinant Human Bak protein ab114338

画像数 1

### 製品の概要

製品名	Recombinant Human Bak protein
タンパク質長	Protein fragment

### 製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	<a href="#">Q16611</a>
生物種	Human
配列	LQPTAENAYEYFTKIATSLFESGINWGRVVALLGFGYRLALHVVYQHGLTG FLGQVTRFVDFMLHHCIARWIAQRGGWVAALNLGNGPILNVLVVLGVVL LGQFVRRFFKS
分子量	38 kDa including tags
領域	100 to 211

### 特性

Our [Abpromise guarantee](#) covers the use of **ab114338** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA SDS-PAGE Western blot
製品の状態	Liquid
備考	Protein concentration is above or equal to 0.05 mg/ml. Best use within three months from the date of receipt of this protein.

### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
------------	---


pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

## 関連情報

機能	In the presence of an appropriate stimulus, accelerates programmed cell death by binding to, and antagonizing the anti-apoptotic action of BCL2 or its adenovirus homolog E1B 19k protein. Low micromolar levels of zinc ions inhibit the promotion of apoptosis.
組織特異性	Expressed in a wide variety of tissues, with highest levels in the heart and skeletal muscle.
配列類似性	Belongs to the Bcl-2 family.
ドメイン	Intact BH3 motif is required by BIK, BID, BAK, BAD and BAX for their pro-apoptotic activity and for their interaction with anti-apoptotic members of the Bcl-2 family.
細胞内局在	Mitochondrion membrane.

## 画像



ab114338 analysed on a 12.5% SDS-PAGE gel stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors