

Product datasheet

Recombinant Human CCL27 protein ab93953

画像数 1

製品の概要

製品名	Recombinant Human CCL27 protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli
アミノ酸配列	
生物種	Human
配列	MFLLPSTAC CTQLYRKPLS DKLLRKVIQV ELQEADGDCH LQAFVLHLAQ RSICIHQNP SLSQWFEHQE RKLHGTLPKL NFGMLRKMG

特性

Our [Abpromise guarantee](#) covers the use of **ab93953** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Mass Spectrometry SDS-PAGE
精製度	> 95 % SDS-PAGE. ab93953 was purified using conventional chromatography techniques.
製品の状態	Liquid
備考	Endotoxin Level: < 1.0 EU per 1 µg of protein (determined by LAL method)


前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Preservative: None Constituents: 10% Glycerol, 10mM Sodium citrate, pH 3.5
------------	---

関連情報

機能	Chemotactic factor that attracts skin-associated memory T-lymphocytes. May play a role in mediating homing of lymphocytes to cutaneous sites. Binds to CCR10.
組織特異性	Testis, thymus, placenta, ovary and skin.
配列類似性	Belongs to the intercrine beta (chemokine CC) family.
細胞内局在	Secreted.

Recombinant Human CCL27 protein 画像



15% SDS-PAGE analysis of ab93953 (3µg)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors