

Product datasheet

Recombinant Human CD62L protein ab53383

画像数 1

製品の概要

製品名	Recombinant Human CD62L protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Escherichia coli

アミノ酸配列

生物種	Human
配列	WT YHYSEKPMNW QRARRFCRDN YTDLVAIQNK AEIEYLEKTL PFSRSYYWIG IRKIGGIWTW VGTNKSLTEE AENWGDGEPN NKKNKEDCVE IYIKRNKDAG KWND DACHKL KAALCYTASC QPWSCSGHGE CVEIINNYTC NCDVGYYPQ CQFVIQCEPL EAPELGTMDC THPLGNFSFS SQCAFSCSEG TNLTGIEETT CGPFGNWSSP EPTCQVIQCE PLSAPDLGIM NCSHPLASFS FTSACTFICS EGTELIGKKK TICESSGIWS NPSPICQKLD KS

領域	39 to 332
----	-----------

特性

Our [Abpromise guarantee](#) covers the use of **ab53383** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
精製度	> 95 % SDS-PAGE.
製品の状態	Liquid

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

Preservative: None

Constituents: 50% Glycerol, 20mM Tris, pH 8

関連情報

機能

Cell surface adhesion protein. Mediates the adherence of lymphocytes to endothelial cells of high endothelial venules in peripheral lymph nodes. Promotes initial tethering and rolling of leukocytes in endothelia.

組織特異性

Expressed in B cell lines and T lymphocytes.

配列類似性

Belongs to the selectin/LECAM family.

Contains 1 C-type lectin domain.


Contains 1 EGF-like domain.

Contains 2 Sushi (CCP/SCR) domains.

細胞内局在

Membrane.

画像



SDS page analysis of ab53383

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors