

Product datasheet

Recombinant Human CD94 protein ab114885

画像数 1

製品の概要

製品名	Recombinant Human CD94 protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Wheat germ

アミノ酸配列

アクセッション番号 [Q13241](#)

生物種 Human

配列
 MAVFKTTLWRLISGTLGIICLSLMATLGILLKNSFTKLSIEPAFTPGPNI
 ELQKSDCCSCQEKWVGYRCNCYFISSEKQTNWESRHLCASQKSSLLQLQ
 NTDELDFMSSSRQFYWIGLSYSEEHTAWLWENGALSQYLFPSFETFNTK
 NCIAYNPNGNALDESCEDKNRYICKQQLI

分子量 45 kDa including tags

領域 1 to 179

特性

Our [Abpromise guarantee](#) covers the use of **ab114885** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	SDS-PAGE
	Western blot

製品の状態 Liquid

備考 Protein concentration is above or equal to 0.05 mg/ml.
 Best use within three months from the date of receipt of this protein.

前処理および保存

保存方法および安定性 Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.


pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

関連情報

機能	Plays a role as a receptor for the recognition of MHC class I HLA-E molecules by NK cells and some cytotoxic T-cells.
組織特異性	Natural killer cells.
配列類似性	Contains 1 C-type lectin domain.
細胞内局在	Membrane.

画像



12.5% SDS-PAGE showing ab114885 at approximately 45.10kDa.
Stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors