

## Product datasheet

# Recombinant Human Cpn10 protein ab91566

### 製品の概要

---

製品名	Recombinant Human Cpn10 protein
タンパク質長	Full length protein

### 製品の詳細

---

由来	Recombinant
由来	Escherichia coli

### アミノ酸配列

生物種	Human
-----	-------

### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab91566** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE
精製度	> 90 % SDS-PAGE. ab91566 was purified using multi-step chromatography
製品の状態	Liquid

### 前処理および保存

---

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Preservative: None Constituents: 20mM Tris, 150mM Sodium chloride, pH 8
------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### 関連情報

---

機能	Eukaryotic CPN10 homolog which is essential for mitochondrial protein biogenesis, together with CPN60. Binds to CPN60 in the presence of Mg-ATP and suppresses the ATPase activity of the latter.
----	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

配列類似性	Belongs to the GroES chaperonin family.
細胞内局在	Mitochondrion matrix.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors