

Product datasheet

Recombinant Human DCUN1D2 protein ab113157

画像数 1

製品の概要

製品名	Recombinant Human DCUN1D2 protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli

アミノ酸配列

アクセッション番号 [Q6PH85](#)

生物種 Human

配列 **MGSSHHHHHSSGLVPRGSHMHKLGKSSQKDKVRQFMACTQAGERTAIYCL**
 TQNEWRLDEATDSFFQNPDSLHRESMRNAVDDKKLERLYGRYKDPQDENK
 IGVDGIQQFCDDLSDPASISVLVIAWKFRAATQCEFSRKEFLDGMTELG
 CDSMEKLLKALLPRLEQELKDTAKFKDFYQFTFTFAKNPGQKGLDLEMAVA
 YWKLVLVSGRFKFLDLWNTFLMEHHKRSIPRDTWNLLLDGFMNIADDMSNY
 DEEGAWPVLIDDFVEYARPVVTGGKRSLF

分子量	32 kDa including tags
領域	1 to 259
タグ	His tag N-Terminus

特性

Our [Abpromise guarantee](#) covers the use of **ab113157** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Mass Spectrometry SDS-PAGE
質量分析	MALDI-TOF
精製度	> 95 % SDS-PAGE. ab113157 was purified using conventional chromatography techniques.
製品の状態	Liquid

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.32% Tris HCl, 0.02% DTT, 30% Glycerol, 0.58% Sodium chloride

関連情報

機能


May contribute to neddylation of cullin components of SCF-type E3 ubiquitin ligase complexes. Neddylation of cullins play an essential role in the regulation of SCF-type complexes activity.

配列類似性

Contains 1 DCUN1 domain.

Contains 1 UBA-like domain.

画像



15% SDS PAGE analysis of 3 µg ab113157.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors