

## Product datasheet

# Recombinant Human EBF1 protein ab112729

### 製品の概要

|        |                                |
|--------|--------------------------------|
| 製品名    | Recombinant Human EBF1 protein |
| タンパク質長 | Protein fragment               |

### 製品の詳細

|    |                  |
|----|------------------|
| 由来 | Recombinant      |
| 由来 | Escherichia coli |

### アミノ酸配列

|           |                        |
|-----------|------------------------|
| アクセッション番号 | <a href="#">Q9UH73</a> |
| 生物種       | Human                  |
| 分子量       | 36 kDa                 |

### 特性

Our [Abpromise guarantee](#) covers the use of **ab112729** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

|          |              |
|----------|--------------|
| アプリケーション | Western blot |
|----------|--------------|

|       |        |
|-------|--------|
| 製品の状態 | Liquid |
|-------|--------|

### 前処理および保存

|            |                                                                                                                                                                                                   |
|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 保存方法および安定性 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.<br>pH: 6.80<br>Constituents: 1.58% Tris HCl, 3.09% DTT, 4% SDS, 0.2% Bromophenol blue, 20% Sucrose |
|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### 関連情報

|    |                                                                                                          |
|----|----------------------------------------------------------------------------------------------------------|
| 機能 | Transcriptional activator which recognizes variations of the palindromic sequence 5'-ATTCCCNNGGGAATT-3'. |
|----|----------------------------------------------------------------------------------------------------------|

|       |                                                          |
|-------|----------------------------------------------------------|
| 配列類似性 | Belongs to the COE family.<br>Contains 1 IPT/TIG domain. |
|-------|----------------------------------------------------------|

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

---

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors