

Product datasheet

Recombinant Human EEF1D protein ab104755

画像数 1

製品の概要

製品名	Recombinant Human EEF1D protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli

アミノ酸配列

アクセッション番号 [P29692](#)

生物種 Human

配列 **MGSSHHHHHSSGLVPRGSHMATNFLAHEKIWFDKFKYDDAERRFYEQMN
GPVAGASRQENGASVILRDIARARENIQKSLAGSSGPGASSGTSGDHGEL
VVRIASLEVENQSLRGVVQELQQAISKLEARLNVLEKSSPGHRATAPQTQ
HVSPMRQVEPPAKKPAEDDEDDIDLFGSDNEEDKEAAQLREERLR
QYAEKKAKKPALVAKSSILLDVKPWDETDMAQLEACVRSIQDGLVWGA
SKLVPVGYGIRKLQIQCVVEDDKVGTDLLEEEITKFEHVQSVDAAFNK
I**

分子量 33 kDa including tags

領域 1 to 281

タグ His tag N-Terminus

特性

Our [Abpromise guarantee](#) covers the use of **ab104755** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Mass Spectrometry

SDS-PAGE

質量分析 MALDI-TOF

精製度 > 90 % SDS-PAGE.

ab104755 is purified using conventional chromatography techniques.

製品の状態 Liquid

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

Preservative: None

Constituents: 10% Glycerol, 20mM Tris HCl, 1mM DTT, pH 8.0

関連情報

機能

Isoform 1: EF-1-beta and EF-1-delta stimulate the exchange of GDP bound to EF-1-alpha to GTP, regenerating EF-1-alpha for another round of transfer of aminoacyl-tRNAs to the ribosome. Isoform 2: Regulates induction of heat-shock-responsive genes through association with heat shock transcription factors and direct DNA-binding at heat shock promoter elements (HSE).

組織特異性

Isoform 2 is specifically expressed in brain, cerebellum and testis.


配列類似性

Belongs to the EF-1-beta/EF-1-delta family.

細胞内局在

Nucleus.

画像



15% SDS-PAGE (3μg).

Note: Molecular weight on SDS-PAGE will appear higher than predicted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors