

Product datasheet

Recombinant Human Glycerol 3 Phosphate Dehydrogenase protein ab93724

画像数 1

製品の概要

製品名	Recombinant Human Glycerol 3 Phosphate Dehydrogenase protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli

アミノ酸配列

生物種	Human
配列	MASKKVCIVG SGNWGSIAIK IVGGNAAQLA QFDPRTMWV FEEDIGGKKL TEIINTQHEN VKYLPGHKLP PNVVAVPDVV QAAEDADILI FVVPHQFIGK ICDQLKGHLK ANATGISLIK GVDEGPNGLK LISEVIGERL GIPMSVLMGA NIASEVADEK FCETTIGCKD PAQGQLLCEL MQTPNFRITV VQEVDTVEIC GALKNVVAVG AGFCDGLGFG DNTKAAVIRL GLMEMIAFAK LFCSGPVSSA TFLESCGVAD LITTCYGGRN RKVAEAFART GKSIEQLEKE LLNGQKLQGP ETARELYSIL QHKGLVDKFP LFMAVYKVCY EQPVGFEFIH CLQNHPEHM

特性

Our [Abpromise guarantee](#) covers the use of **ab93724** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
精製度	> 90 % SDS-PAGE. ab93724 is purified using conventional chromatography techniques.
製品の状態	Liquid

## 前処理および保存

### 保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

Preservative: None

Constituents: 10% Glycerol, 20mM Tris HCl, pH 8.0

## 関連情報


### 配列類似性

Belongs to the NAD-dependent glycerol-3-phosphate dehydrogenase family.

### 細胞内局在

Cytoplasm.

## 画像



15% SDS-PAGE showing ab93724 at approximately 37.5kDa (3μg).

SDS-PAGE - Recombinant Human Glycerol 3  
Phosphate Dehydrogenase protein (ab93724)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors