

## Product datasheet

# Recombinant Human HDDC3 protein ab113583

### 画像数 1

#### 製品の概要

製品名	Recombinant Human HDDC3 protein
タンパク質長	Full length protein

#### 製品の詳細

由来	Recombinant
由来	Escherichia coli

#### アミノ酸配列

アクセッション番号 [Q8N4P3](#)

生物種 Human

配列 **MGSSHHHHHSSGLVPRGSHLDEVELHFGAQVRRLVEEVTDDKTLPKLER**  
**KRLQVEQAPHSSPGAKLVKLADKLYNLRDLNRCTPEVKIQMGSEAAQLLE**  
**AADFAARKHRQRRKDPEGTPYINHPIGVARILTHEAGITDIVLQAALL**  
**HDTVEDTDTT**

分子量 18 kDa including tags

領域 1 to 140

タグ His tag N-Terminus

#### 特性

Our [Abpromise guarantee](#) covers the use of **ab113583** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
精製度	> 90 % SDS-PAGE. ab113583 is purified using conventional chromatography techniques: anion exchange chromatography followed by gel filtration chromatography with 20mM Tris pH 7.5, 2mM EDTA. No detergents such as urea, Triton or Tween were used to purify the proteins. They purified the proteins by using
製品の状態	Liquid

## 前処理および保存

### 保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.32% Tris HCl, 40% Glycerol, 0.88% Sodium chloride, 0.02% DTT

## 関連情報

### 機能


ppGpp hydrolyzing enzyme involved in starvation response.

### 配列類似性

Belongs to the MESH1 family.

Contains 1 HD domain.

## 画像



15% SDS-PAGE showing ab113583 at approximately 17.9 kDa (3µg).

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors