

## Product datasheet

# Recombinant Human HDHD3 protein ab125485

### 製品の概要

製品名	Recombinant Human HDHD3 protein
タンパク質長	Full length protein

### 製品の詳細

由来	Recombinant
由来	Escherichia coli

### アミノ酸配列

アクセッション番号 [Q9BSH5](#)

生物種 Human

配列  
 MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD  
 KDRWGSMHR LQIRLLTWDV KDTLLRLRHP  
 LGEAYATKAR AHGLEVEPSA LEQGFRQAYR  
 AQSHSFPNYG LSHGLTSRQW WLDVVLQTFH  
 LAGVQDAQAV APIAEQLYKD FSHPCTWQVL  
 DGAEDTLREC RTRGLRLAVI SNFDRRLEGI  
 LEGLGLREHF DFVLTSEAAG WPKPDPRIFQ  
 EALRLAHMEP VVAHVGDNY LCDYQGPRAV  
 GMHSFLVVG P QALDPVVRDS VPKEHILPSL  
 AHLLPALDCL EGSTPGL

分子量 32 kDa including tags

領域 1 to 251

タグ His tag N-Terminus

配列の追加情報 HDHD3 Human Recombinant fused with a36 amino acid His tag at N-terminus produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 287 amino acids (1-251 a.a.).

### 特性

Our [Abpromise guarantee](#) covers the use of **ab125485** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

精製度 > 95 % SDS-PAGE.

ab125485 is purified by proprietary chromatographic techniques.

#### 製品の状態

Liquid

#### 備考

For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

#### 前処理および保存

---

##### 保存方法および安定性

Shipped at 4°C. Store at -20°C. Store under desiccating conditions.

pH: 8.00

Constituents: 0.32% Tris HCl, 10% Glycerol

#### 関連情報

---

##### 配列類似性

Belongs to the HAD-like hydrolase superfamily.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors