

Product datasheet

Recombinant Human LHPP protein ab116175

画像数 1

製品の概要

製品名	Recombinant Human LHPP protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli

アミノ酸配列

アクセッション番号 [Q9H008](#)

生物種 Human

配列 MRGSHHHHP WYASMTGGQQ MGRDLYDDDD  
 KDRWGSHMAP WGKRLAGVRG VLLDISGVLY  
 DSGAGGGTAI AGSVEAVARL KRSRLKVRFC  
 TNESQKSRAE LVGQLQRLGF DISEQEV TAP  
 APAACQILKE QGLRPYLLIH DGVRSEFDQI  
 DTSNPNCVVI ADAGESFSYQ NMNNAFQVLM  
 ELEKPV L ISL GKGRYYKETS GLMLDVGPYM  
 KALEYACGIK AEVVGKPSPE FFKSALQAIG  
 VEAHQAVMIG DDIVGDVGG A QRCGMRALQV  
 RTGKFRPSDE HHPEVKADGY VDNLAEAVDL LLQHADK

分子量 34 kDa including tags

領域 1 to 270

タグ His tag N-Terminus , DDDDK tag N-Terminus

特性

Our [Abpromise guarantee](#) covers the use of **ab116175** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Mass Spectrometry
質量分析	MALDI-TOF

精製度	> 90 % SDS-PAGE. ab116175 was purified using conventional chromatography.
製品の状態	Liquid


#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  pH: 8.00 Constituents: 10% Glycerol, 0.32% Tris HCl, 0.02% DTT, 0.58% Sodium chloride
------------	--

#### 関連情報

組織特異性	Expressed in liver, kidney and moderately in brain.
配列類似性	Belongs to the HAD-like hydrolase superfamily.

#### 画像



ab116175 on a 15% SDS-PAGE (3ug)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors